

BLOTTING PAPERS

Recommended for blotting techniques

- Made from ultrapure chromatography paper
- High and uniform absorbency

Technical specifications

Grade	Properties	Surface	Thickness mm	Applications
BP3469	Medium Absorbency	Smooth	0.4	Wicking & Gel support
BP2668	Medium Absorbency	Smooth	0.8	Membrane Gel support
BP2294	High Absorbency	Smooth	1.2	Semidry Blotting

(1) Measures in 30 min.

* measures in 10 min

Applications

BP3469


Southern, Northern and Western blots
Dot and slot-blots
Lifting of sequencing else
Buffer wicking
General laboratory use

BP2668

Lysis / denaturation of colony or plaque lifts
Western blots

BP2294

Semidry blotting of proteins



Ordering information

Formats and size on request. Please, contact us