

## **FILTER PAPERS FOR THE SUGAR INDUSTRY**

**Recommended for the clarifying filtration of beet extracts and juices**

- High filtration speed combined with good retention of particles
- 2 surfaces available: smooth and creped

### **Technical Specifications**

<b>Grade</b>	<b>Properties</b>	<b>Filtration Herzbergs</b>	<b>Weight g/m<sup>2</sup></b>	<b>Thickness mm</b>
<b>FP3000</b>	Fast, smooth	95	69	0.16
<b>FP3002</b>	Medium fast , smooth	150	61	0.14
<b>FP3003</b>	Medium Fast, creped	240	65	0.25
<b>FP3459</b>	Fast, creped	110	74	0.30
<b>FP3533</b>	Fast, creped	80	72	0.27

### **Applications**

- Clarification of dried beet pulp extracts
- Filtration of beet juice after addition of lead acetate for polarimetric sugar determination
- FP3459 and FP3533 recommended for Venema unit according with the lead acetate method.

### **Ordering information**

<b>Grade</b>	<b>FP3000</b>	<b>FP3002</b>	<b>FP3003</b>	<b>FP3459</b>	<b>FP3533</b>
<b>Size</b>					
<b>185 mm</b>	DP3000185	DP3002185	DP3003185		
<b>200 mm</b>	DP3000200	DP3002200	DP3003200		
<b>230 mm</b>				DP3459230	DP3533230
<b>250 mm</b>	DP3000250	DP3000250	DP3003250		

DP: Round filters

Presentation : 1000 units/box

**Other size available on request**

